

Instruction for Use

【Product Name】

Gauze Swabs

【Product Composition】

The product is mainly composed of white or dyed bleached gauze, gauze is made of 100% cotton.

【Intended Use】

It's intended to be used as mechanical barrier, for compression or for absorption of exudates. This is a single use device.

【Intended User】

Medical staff or lay person

【Intended patient population】

Used for all age groups except infants.

【Handling Instruction】

- Clean the wound area, hands and instruments in accordance with normal procedures.
- Check the package, if the package is not intact, or the product is abnormal compared with others, it must be replaced.
- Open the packaging and take out the product, gently place the product on the wound that needs to absorb exudates.

【Warnings and Precautions】

- The product is disposable and cannot be reused.
- Stop using if product has any defect.
- Stop using if the product has expired.

Document Name: Instruction for Use

Document No.: IFU-GS-M020105

Product : Gauze Swabs

Version: 01



【Shelf Life】 5years

【Contraindications】

Do not use on patients with a known hypersensitivity to the product itself or to its components.

【Side Effects】

None

【Disposal】

Disposal should be handled according to local environment procedures.

【Special Storage Conditions and Methods】



















The product should be stored in a cool and dry place and avoid direct sunlight.

The transportation should prevent the weight, sun and rain.

【Notice】

Any serious adverse event related to this product must be report to Allmed Medical Products Co., Ltd or authorized representative or national health authority in your country.

【Explanation of figures, symbols and abbreviations used in labels】

| | | | | | | | |
|---|----------------------------------|---|-------------------------|--|---|--|------------------------------|
|  | Do not reuse |  | EO sterilization |  | Warnings and Precautions |  | Manufacturer Name Address |
|  | Used-by date |  | CE Symbol |  | Manufacture Date |  | Consult instructions for use |
|  | Batch Code |  | Catalogue Number |  | Name and Address of European Union Representative |  | medical device |
|  | Single sterile barrier system |  | Keep away from sunlight |  | Keep dry |  | Do not resterilize |
|  | Do not use if package is damaged |  | importer | | | | |

【Manufacturer】 Allmed Medical Products Co.,Ltd

【Address】 No.18 Qi xing Road, Majiadian Town, 443200 Zhijiang City, Hubei Province , CHINA

【European Representative】 Shanghai International Holding Corp. GmbH(Europe)

【Address】 Eiffestrasse 80, 20537 Hamburg, Germany

【Revision Date】 2024-09-20

Document Name: Instruction for Use

Document No.: IFU-GS-M020105

Product : Gauze Swabs

Version: 01



ALLMED MEDICAL PRODUCTS CO.LTD