


175x115x175mm

Die line(刀线) : 
Pressure line(压线) : 
perforated line(针孔刀) : 

5mm :2mm
5mm :2mm
5mm :2mm



| Body Part | Gauge |
|---|-------|
| Fingers (hand and foot), wrist and ankle | 0-1 |
| Ankle, foot, elbow, arm and forearm | 2 |
| Arm, forearm, leg and knee | 3 |
| Shoulder, head and thigh (pediatric use) | 4 |
| Head, thigh and armpit (infant use) | 5 |
| Perineum, armpit, pelvis, thorax, trunk, head and thigh | 6 |

Refer to the leaflet for guidance on use and bandage applications.

General Observations

Healing material, not sterilized. Store in a cool, dry place. Do not expose to sunlight.

Manufactured with Polyamide (70%) and Latex Rubber (30%). The tubular net may be applied depending on the wound conditions and if the product allows it.

Attention: it is important to evaluate the wound's needs; the tubular net is recommended for use on superficial wounds.



Elastic mesh bandage. Tubular shape.

100 mts

Caliber. #5

Lot no.:
Manufacturing Date:
Expiration Date:



Made in China: ALLMED MEDICAL PRODUCTS CO., LTD.
No. 16 Qixing Road, Majadian Town, Zhijiang City, Hubei Province, 443200, China.