



**Transparent, microporous,  
self-adhesive dressing**

**Measures:  
7.0 to 8.5 cm × 5.08 to 6.0 cm**

**Instructions for use:**

1. Open the package and remove the dressing.
2. Peel back the protective film to expose the adhesive surface.
3. Place the dressing over the site and press to secure it.
4. Remove the front frame or protective film while smoothing the edges.
5. Complete and apply the label if included.

**Lot no.:**  
**Manufacturing Date:**  
**Expiration Date:**

 **Sterile  
and disposable**

\*The sterility of the dressing cannot be guaranteed if the sachet shows signs of having been previously ruptured.  
\*Store in a cool, dry place.

**50 pieces**

**STERILE** **EO** **CE** 

**UDI** 6941558073174

**Made in China: ALLMED MEDICAL PRODUCTS CO., LTD.**  
No. 18 Qixing Road, Majiadian Town, Zhijiang City, Hubei Province, 443200, China.