

330x65x330mm

Die line(刀线) :
Pressure line(压线) :
perforated line(冲孔刀) :



Body Part	Gauge
Fingers (hand and foot), wrist and ankle	0-1
Ankle, foot, elbow, arm and forearm	2
Arm, forearm, leg and knee	3
Shoulder, head and thigh (pediatric use)	4
Head, thigh and armpit (infant use)	5
Perineum, armpit, pelvis, thorax, trunk, head and thigh	6


Refer to the leaflet fot guidance on use and bandage applications.

General Observations

Healing material, not sterilized. Store in a cool, dry place. Do not expose to sunlight.

Manufactured with Polyamide (70%) and latex Rubber (30%). The tubular net may be applied depending on the wound conditions and if the product allows it.

Attention: it is important to evaluate the wound´s needs; the tubular net is recommended for use on superticial wounds.



Elastic mesh bandage.
Tubular shape.
100 mts
Caliber. #6

Lot no:
Manufacturing
Date:
Expiration
Date:



CE UDI 6941558073082NSZ6

Made in China: ALLMED MEDICAL PRODUCTS CO., LTD.
No.18 Qixing Road, Majiadian Town, Zhijiang City, Hubei Province, 443200, China.