



## Dressing

**Sterile, transparent polyurethane film with a hypoallergenic, latex-free adhesive and Comfort Technology; the film includes a transparent, absorbent gel pad impregnated with 2% chlorhexidine gluconate.**

**7.0 x 8.5 cm. Pad: 2.0 x 2.0 cm**

Lot no.:  
Manufacturing Date:  
Expiration Date:

 **Sterile  
and disposable**

The sterility of the dressing cannot be guaranteed if the sachet shows signs of having been previously ruptured.  
Store in a cool, dry place.

**1 piece**

**STERILE** **E0**  

**Made in China:** ALLMED MEDICAL PRODUCTS CO., LTD.  
No. 18 Qixing Road, Majiadian Town, Zhijiang City, Hubei Province, 443200, China.